



Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





MYC/IgH/Cen 8 t(8;14)

Methodology

FISH

Test Description

Probes: Trisomy 8 (Cen 8) | MYC/IgH t(8;14)

Disease(s): Burkitt lymphoma, NHL

Clinical Significance

The MYC/IgH rearrangement is common in Burkitt lymphoma and diffuse large B-cell lymphoma, associated with a more aggressive tumor phenotype and a lower survival rate. MYC/IgH/CEN8 can confirm heavy chain rearrangement when MYC is rearranged.

Specimen Requirements

- **Bone marrow aspirate:** 1-2 mL in sodium heparin tube. EDTA tube is acceptable.
- **Peripheral blood:** 2-5 mL in sodium heparin tube. EDTA tube is acceptable.
- **Fresh tissue:** Tissue in RPMI.
- **Bone Marrow/ Peripheral Blood Smear or Fresh Tissue Touch Preparation Slides:** minimum 1 slide labeled with specimen type.
- **Fluids:** Equal parts RPMI to specimen volume.
- **Fixed Cell Suspension:** A client fixed cell suspension may be submitted for testing as long as it is received in 3:1 Methanol:Glacial Acetic Acid.
- **Paraffin block:** H&E slide (required) plus paraffin block. Circle H&E for tech-only.
- **Cut slides:** H&E slide (required) plus 2 unstained slides cut at 4 microns. Circle H&E for tech-only.
- **Note:** Please exclude biopsy needles, blades, and other foreign objects from transport tubes. These can compromise specimen viability and yield, and create hazards for employees.

Storage & Transportation

Refrigerate specimen. Do not freeze. Use cold pack for transport, making sure cold pack is not in direct contact with specimen. For fresh samples: ship same day as drawn whenever possible; specimens <72 hours old preferred.

CPT Code(s)*

88374x1 automated. Codes may differ if manual analysis is performed.

New York Approved

Yes

Level of Service

Technical, Global

Turnaround Time

3-5 days for unfixed or FFPE specimens

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way
Fort Myers, FL 33912
Phone: 239.768.0600/ Fax: 239.690.4237
neogenomics.com
© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved.
All other trademarks are the property of their respective owners
Rev. 042924